



Product Datasheet

Lysozyme antibody (orb334039)


Description	Mouse monoclonal antibody to Lysozyme
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	CyTOF, FC, ICC
Immunogen	human lysozyme
Target	Lysozyme
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Intracellular staining; recommended dilution: 1-12 µg/ml. Excellent. Immunohistochemistry (paraffin sections): Positive tissue: tonsil.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P61626
Entrez	4069
Dilution Range	Flow cytometry: Intracellular staining; recommended dilution: 1-12 µg/ml. Excellent. Immunohistochemistry (paraffin sections): Positive tissue: tonsil.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (intracellular staining) S...



Mass cytometry (intracellular staining) ...



Immunocytochemistry staining of lysozyme...