


Product Datasheet

PRKACA antibody (orb334051)

Description	Goat polyclonal antibody to PRKACA
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRKACA
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 48kDa band observed in lysates of cell line HeLa and MCF7 (calculated MW of 49.1kDa according to NP_002721.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	ESPAQNTAHL DQ
Clonality	Polyclonal
MW	49.1; 40.6; 39.8
Entrez	5566
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa (Lane 1) a...