

Product Datasheet

NLRP3 antibody (orb334052)

Description	Goat polyclonal antibody to NLRP3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Cryopirin / NALP3 / NLRP3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Bone Marrow and Peripheral Blood Lymphocytes lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been re
Protein Sequence	NQERTSYLEKKLSCK
Clonality	Polyclonal
MW	118; 106; 112; 112
Entrez	114548
Dilution Range	ELISA: 1:4000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.