


Product Datasheet

UBE2I antibody (Biotin) (orb334070)

Description	Goat polyclonal antibody against UBE2I conjugated to Biotin.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Target	UBE2I, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 17kDa band observed in Human Kidney lysates (calculated MW of 18.0kDa according to Human NP_003336.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.3-1µg/ml.. IHC: .
Protein Sequence	SGIALSRLAQERKC.
Clonality	Polyclonal
MW	18.0; 20.5
Entrez	7329
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...