

Product Datasheet

MYD88 antibody (Biotin) (orb334076)

Description Goat polyclonal antibody against MYD88 conjugated to Biotin.

Species/Host Goat

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Target MYD88, Biotinylated

Preservatives Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage Aliquot and store at -20°C. Minimize freezing and thawing.

Note For research use only

Application notes WB: Approx 33kDa band observed in Human Thymus lysates (calculated MW of 33.2kDa according to NP_002459.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.1-0.3µg/ml.

Protein Sequence IKYKAMKKEFP.


Clonality Polyclonal

MW 35.4; 34.6; 29.6; 22.1; 17.2; 28.5; 23.5

Entrez [4615](#)

Dilution Range ELISA: 1:16000, WB: 0.1-0.3 µg/ml

Expiration Date 12 months from date of receipt.



Western blot analysis of Human Thymus ly...