


Product Datasheet

MS4A1 antibody (Biotin) (orb334081)

Description	Goat polyclonal antibody against MS4A1 conjugated to Biotin.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CD20 / MS4A1 (C Terminus), biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 35kDa band observed in lysates of cell line Daudi (calculated MW of 33.1kDa according to NP_068769.2). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.3-1µg/ml.
Protein Sequence	CQDQESSPIENDSSP.
Clonality	Polyclonal
MW	33.1
Entrez	931
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Daudi lysate us...