


Product Datasheet

PVALB antibody (Biotin) (orb334085)

Description	Goat polyclonal antibody against PVALB conjugated to Biotin.
Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Parvalbumin, biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 12kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 12.1kDa according to NP_002845.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 1-3µg/ml.
Protein Sequence	GVDEFSTLVAES.
Clonality	Polyclonal
MW	12.1; 11.7
Entrez	5816
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...