



## Product Datasheet

### Transthyretin antibody (orb334220)


|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Transthyretin antibody                                                                                         |
| <b>Species/Host</b>        | Mouse                                                                                                          |
| <b>Reactivity</b>          | Human, Rat                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                   |
| <b>Tested Applications</b> | ELISA, IHC-P, WB                                                                                               |
| <b>Immunogen</b>           | Human Transthyretin recombinant protein (orb315545) was used as immunogen.                                     |
| <b>Target</b>              | Transthyretin. Recognizes endogenous levels of Transthyretin protein                                           |
| <b>Form/Appearance</b>     | Liquid: 10 mM PBS, 0.02% sodium azide                                                                          |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml                                                 |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                          |
| <b>Isotype</b>             | IgG1                                                                                                           |
| <b>Clonality</b>           | Monoclonal                                                                                                     |
| <b>Purity</b>              | Monoclonal antibodies are purified by peptide affinity chromatography                                          |
| <b>MW</b>                  | 16 kDa                                                                                                         |
| <b>Uniprot ID</b>          | <b>P02767, P02766, P07309</b>                                                                                  |
| <b>NCBI</b>                | <b>000362, 000371</b>                                                                                          |
| <b>Entrez</b>              | <b>7276</b>                                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                |



Western blot analysis of His-28TTR recom...



Western blot analysis of His-28TTR recom...



IHC-P image of rat pancreas tissue using...