



Product Datasheet

alpha 1 Antitrypsin antibody (orb334222)


Description	alpha 1 Antitrypsin antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide (please contact us for details)
Target	alpha 1 Antitrypsin. Recognizes endogenous levels of alpha 1 Antitrypsin protein
Form/Appearance	Liquid: 10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Monoclonal antibodies are purified by peptide affinity chromatography
MW	46 kDa
Uniprot ID	P17475, Q00896, P01009
NCBI	000295, 000286
Entrez	5265
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of Trx - 32A1AT re...



IHC-P image of rat stomach tissue using
...