

Product Datasheet

Fish auratus Leptin protein (orb334231)

Description	Recombinant Fish auratus Leptin protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	8His tag at the N-terminus.
Protein Sequence	Pro22-Cys171
Source	P.Pichia
Uniprot ID	B8YI02
Expiration Date	6 months from date of receipt.