

Product Datasheet

E. coli Tryptophan Synthase alpha Chain protein (orb334237)

Description	Recombinant E. coli Tryptophan Synthase alpha Chain protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Ser268
Source	E.coli
Uniprot ID	P0A877
Expiration Date	6 months from date of receipt.