

Product Datasheet

Human EpCAM protein (orb334261)

Description	Recombinant human EpCAM protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-Fc
Note	For research use only
Application notes	Fc tag at the C-terminus.
Protein Sequence	Gln24-Lys265
Source	Human Cells
Uniprot ID	AAH14785.1
Expiration Date	6 months from date of receipt.