

Product Datasheet

Human Isocitrate Dehydrogenase 1 protein (orb334293)

Description	Recombinant human Isocitrate Dehydrogenase 1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM Tris, 200mM NaCl, 10% glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM Tris-HCl, 200mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	8His tag at the C-terminus.
Protein Sequence	Met1-Leu414
Source	E.coli
Uniprot ID	O75874
Expiration Date	6 months from date of receipt.