

Product Datasheet

Human beta-1,3-Glucuronyltransferase 3 protein (orb334334)

Description	Recombinant human beta-1,3-Glucuronyltransferase 3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaC,2mM EDTA,20% Glycerol,pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaC, 2mM EDTA, 20% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Glu72-Val335
Source	E.coli
Uniprot ID	O94766
Expiration Date	6 months from date of receipt.