

Product Datasheet

Mouse Fms-Like Tyrosine Kinase 3 Ligand protein (orb334374)

Description	Recombinant mouse Fms-Like Tyrosine Kinase 3 Ligand protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Application notes	Fc tag at the C-terminus.
Protein Sequence	Gly27-Arg188
Source	Human Cells
Uniprot ID	P49772
Expiration Date	6 months from date of receipt.