

Product Datasheet

Mouse Intercellular Adhesion Molecule 2 protein (orb334395)

Description	Recombinant mouse Intercellular Adhesion Molecule 2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Application notes	Fc tag at the C-terminus.
Protein Sequence	Lys23-Gln222
Source	Human Cells
Uniprot ID	P35330
Expiration Date	6 months from date of receipt.