

Product Datasheet

Mouse Interleukin-18 protein (orb334406)

Description	Recombinant mouse Interleukin-18 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mMTris-HCl,1mMEDTA,2mMDTT,5% Trehalose,pH8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Asn36-Ser192
Source	E.coli
Uniprot ID	P70380
Expiration Date	6 months from date of receipt.