

Product Datasheet

Mouse Interleukin-22 protein (orb334410)

Description	Recombinant mouse Interleukin-22 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB,150mM NaCl,, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-mFc
Note	For research use only
Application notes	mFc tag at the C-terminus.
Protein Sequence	Leu34-Val179
Source	Human Cells
Uniprot ID	Q9JJY9
Expiration Date	6 months from date of receipt.