

Product Datasheet

Mouse Interleukin-4 protein (orb334412)

Description	Recombinant mouse Interleukin-4 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 300mM NaCl, 5% Trehalose, pH 6.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	His23-Ser140
Source	E.coli
Uniprot ID	P07750
Expiration Date	6 months from date of receipt.