

Product Datasheet

Mouse Serpin F2 protein (orb334441)

Description	Recombinant mouse Serpin F2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH 8.0(20% glycerol,1mM DTT).
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Val28-Lys491
Source	Human Cells
Uniprot ID	Q61247
Expiration Date	6 months from date of receipt.