

## Product Datasheet

# Mouse Tissue Inhibitors of Metalloproteinases 1 protein (orb334451)

<b>Description</b>	Recombinant mouse Tissue Inhibitors of Metalloproteinases 1 protein
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	tag free
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Cys25-Arg205
<b>Source</b>	Human Cells
<b>Uniprot ID</b>	<b>P12032</b>
<b>Expiration Date</b>	6 months from date of receipt.