

Product Datasheet

Mouse Transforming Growth Factor beta-1 protein (orb334453)

Description	Recombinant mouse Transforming Growth Factor beta-1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50mM Glycine,150mM NaCl,pH 2.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 4mM HCl.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Ala279-Ser390
Source	Human Cells
Uniprot ID	P04202
Expiration Date	6 months from date of receipt.