

Product Datasheet

Mouse alpha-2-HS-Glycoprotein protein (orb334461)

Description	Recombinant mouse alpha-2-HS-Glycoprotein protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	10His tag at the N-terminus.
Protein Sequence	Ala19-Ile345
Source	Human Cells
Uniprot ID	P29699
Expiration Date	6 months from date of receipt.