

Product Datasheet

Vibrio cholerae Toxin B protein (orb334470)

Description	Recombinant Vibrio cholerae Toxin B protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50mM Tris, 200uM NaCl,pH8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM Tris-HCl, 200mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Thr22-Asn124
Source	E.coli
Uniprot ID	P01556
Expiration Date	6 months from date of receipt.