

Product Datasheet

CD153/TNFSF8 Antibody (orb334482)

Description

CD153/TNFSF8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC

Applications
Immunogen

E.coli-derived human CD153 recombinant protein (Position: Q63-D234). Human CD153 shares 75% amino acid (aa) sequence identity with mouse CD153.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

26017 MW

Uniprot ID


P32971

Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, By HeatFlow Cytometry, 1-3µg/1x10⁶cells

Expiration Date

12 months from date of receipt.



CD153 was detected in paraffin-embedded ...



Flow Cytometry analysis of U937 cells us...