


Product Datasheet

TRAP1 Antibody (orb334566)

Description	TRAP1 Antibody	130KD- 1 2 3 4 5 6
Species/Host	Rabbit	70KD-
Reactivity	Human, Rat	55KD-
Conjugation	Unconjugated	35KD-
Tested	IHC, WB	25KD-
Applications		15KD-
Immunogen	E.coli-derived human TRAP1 recombinant protein (Position: A571-H704). Human TRAP1 shares 91.7% and 94% amino acid (aa) sequence identity with mouse and rat TRAP1, respectively.	WB analysis of TRAP1 using anti-TRAP1 an...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	IHC analysis of TRAP1 using anti-TRAP1 a...
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	80110 MW	
Uniprot ID	Q12931	
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Western blot, 0.1-0.5µg/ml, Human, Rat	
Expiration Date	12 months from date of receipt.	