


Product Datasheet

Ceruloplasmin/CP Antibody (orb334586)


Description

Ceruloplasmin/CP Antibody


Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	E. coli-derived mouse Ceruloplasmin recombinant protein (Position: R20-M258). Mouse Ceruloplasmin shares 80.8% and 91.2% amino acid (aa) sequence identity with human and rat Ceruloplasmin, respectively.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	121151 MW
Uniprot ID	Q61147
Dilution Range	Western blot, 0.1-0.5µg/ml, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat, By Heat
Expiration Date	12 months from date of receipt.



WB analysis of Ceruloplasmin/CP using an...



IHC analysis of Ceruloplasmin/CP using a...



IHC analysis of Ceruloplasmin/CP using a...