


Product Datasheet

FCGRT Antibody (orb334600)

Description	FCGRT Antibody	130KD — 100KD — 70KD — 55KD — 35KD — 25KD — 15KD —
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	E. coli-derived mouse FCGRT recombinant protein (Position: Q54-D291). Mouse FCGRT shares 70.2% and 92.4% amino acid (aa) sequence identity with human and rat FCGRT, respectively.	WB analysis of FCGRT using anti-FCGRT an...
Form/Appearance		
Concentration		
Storage		
Note	For research use only	IHC analysis of FCGRT using anti-FCGRT a...
Application notes		
Isotype		
Clonality	Polyclonal	
MW	40093 MW	
Uniprot ID	Q61559	
Dilution Range		
Expiration Date		