

Product Datasheet

NAP2/PPBP Antibody (orb334611)

Description

NAP2/PPBP Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E. coli-derived humn NAP2 recombinant protein (Position: S55-D128).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

13894 MW

Uniprot ID


P02775

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, - Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat ELISA(Cap) , 1-5µg/ml, Human, -

Expiration Date


12 months from date of receipt.



WB analysis of NAP2 using anti-NAP2 anti...



Sandwich ELISA - Recombinant human CXCL7...



IHC analysis of NAP2 using anti-NAP2 ant...