



Product Datasheet

SMYD3 Antibody (orb334626)

Description		
SMYD3 Antibody		
Species/Host	Rabbit	 WB analysis of SMYD3 using anti-SMYD3 an...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC, WB	
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SMYD3 (388-428aa QAMKNLRLAFDIMRVTHGREHSLIEDLILLLEECDANIRAS), different from the related mouse sequence by one amino acid.	
Form/Appearance	Lyophilized	 IHC analysis of SMYD3 using anti-SMYD3 a...
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	49097 MW	
Uniprot ID	Q9H7B4	
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Western blot, 0.1-0.5µg/ml, Human	
Expiration Date	12 months from date of receipt.	