

Product Datasheet

TCF7L1 Antibody (orb334630)

Description

TCF7L1 Antibody

97KD —

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

WB analysis of
TCF7L1 using
anti-TCF7L1 ...

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human TCF7L1 (561-588aa SFPATLHAHQALPVLQAQPLSLVTKSAH), different from the related mouse sequence by one amino acid.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results.
By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

62631 MW

Uniprot ID



[Q9HCS4](#)

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



IHC analysis
of TCF7L1
using anti-
TCF7L1...