

Product Datasheet

NIRF/UHRF2 Antibody (orb334639)

Description

NIRF/UHRF2 Antibody

130KD — 1 2

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human NIRF (15-54aa TIEDVSRKATIEELRERVWALFDVVRPECQRLFYRGKQLEN), identical to the related mouse and rat sequences.

WB analysis of NIRF using anti-NIRF anti...

Form/Appearance


Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.




IHC analysis of NIRF using anti-NIRF ant...

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.



IHC analysis of NIRF using anti-NIRF ant...

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

89985 MW

Uniprot ID
Q96PU4
Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat, By Heat Western blot, 0.1-0.5µg/ml, Human, Rat

Expiration Date

12 months from date of receipt.