


Product Datasheet

Glucocorticoid Receptor antibody (orb334648)

Description	Rabbit polyclonal antibody to Glucocorticoid Receptor
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein corresponding to human Glucocorticoid Receptor.
Target	NR3C1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P06536 , P04150 , P06537
Entrez	24413 , 2908
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of Glucocorticoid

...