



Product Datasheet

HIF1 alpha antibody (orb334653)

Description	Rabbit polyclonal antibody to HIF1 alpha
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant protein corresponding to mouse HIF1 alpha.
Target	HIF1A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q61221 , O35800
Entrez	15251 , 29560
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HIF1 alpha expr...



Immunohistochemical analysis of HIF1 alp...