

Product Datasheet

Histone H1 (MonoMethyl K25) antibody (orb334656)

Description

Rabbit polyclonal antibody to Histone H1

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

Applications

IH, WB

Immunogen

KLH-conjugated synthetic Mono-methylated peptide corresponding to residues surrounding K25 of human Histone H1 protein. The exact sequence is proprietary.

Target

HIST1H1A; HIST1H1C; HIST1H1D; HIST1H1E

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P43275](#), [P16403](#), [Q02539](#), [P10412](#), [P15864](#), [D4A3K5](#), [P43277](#), [P43274](#), [P15865](#), [P16402](#)

Entrez


[291145](#), [3007](#), [3006](#), [3008](#), [50708](#), [80838](#), [50709](#), [3024](#), [14957](#), [201097](#)

Dilution Range


WB: 1:500-1:1000, IHC-P: 1:100-1:200

Expiration Date

12 months from date of receipt.



Western blot analysis of Histone H1 (Monoclonal)



Immunohistochemical analysis of Histone H1 (Monoclonal)