


Product Datasheet

Histone H1 antibody (orb334661)

Description	Rabbit polyclonal antibody to Histone H1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Histone H1. The exact sequence is proprietary.
Target	HIST1H1A; HIST1H1C; HIST1H1D; HIST1H1E
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P43275, P16403, Q02539, P10412, P15864, P43277, P43274, P16402
Entrez	3007, 3006, 3008, 50708, 80838, 50709, 3024, 14957
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H1 expr...



Immunohistochemical analysis of Histone ...