


Product Datasheet

Histone H2A (AcK5) antibody (orb334663)

Description

Rabbit polyclonal antibody to Histone H2A



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IH, WB

Applications


Western blot analysis
of Histone H2A (Ac...

Immunogen

KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K5 of human Histone H2A protein. The exact sequence is proprietary.

Target

HIST1H2AG; HIST1H2AI; HIST1H2AK;
HIST1H2AL; HIST1H2AM; HIST2H2AA3;
HIST2H2AA4; HIST3H2A



Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Immunohistochemical
analysis of Histone ...

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P0CC09](#), [Q6GSS7](#), [P02262](#), [Q7L7L0](#),
[Q8BFU2](#), [Q4FZT6](#), [Q6FI13](#), [P22752](#), [P0C0S8](#)

Entrez

[723790](#), [365877](#), [64646](#), [8329](#), [8330](#), [8332](#),
[319164](#), [8336](#), [15267](#), [92815](#), [8969](#), [319162](#),
[8337](#)

Dilution Range

WB: 1:1000-1:2000, IHC-P: 1:200-1:500

Expiration Date

12 months from date of receipt.