


Product Datasheet

Histone H2A.Z (AcK7) antibody (orb334672)

Description	Rabbit polyclonal antibody to Histone H2A.Z (Acetyl-K7)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K7 of human Histone H2A.Z protein. The exact sequence is proprietary.
Target	H2AFZ
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P0C0S6, P0C0S7, P0C0S5
Entrez	58940, 3015, 51788
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2A.Z (...)



Immunohistochemical analysis of Histone ...