


Product Datasheet

Histone H2B (AcK12) antibody (orb334674)

Description	Rabbit polyclonal antibody to Histone H2B (Acetyl-K12)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K12 of human Histone H2B protein. The exact sequence is proprietary.
Target	HIST1H2BC
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P62807, Q96A08, P70696, Q6ZWY9, P33778, Q00729, Q64475
Entrez	8344, 8346, 24829, 319178, 319179, 8339, 255626, 8347, 319177, 3017, 3018, 8343
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2B (Ac...



Immunohistochemical analysis of Histone ...