

Product Datasheet


Histone H2B (Phospho-S32) antibody (orb334683)

Description

Rabbit polyclonal antibody to Histone H2B



Western blot analysis of Histone H2B (Ph...



Immunohistochemical analysis of Histone ...

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S32 of human Histone H2B protein. The exact sequence is proprietary.
Target	HIST1H2BC
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P62807 , Q96A08 , P70696 , Q6ZWY9 , P33778 , Q00729 , Q64475
Entrez	8344 , 8346 , 24829 , 319178 , 319179 , 8339 , 255626 , 8347 , 319177 , 3017 , 3018 , 8343
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.