


Product Datasheet

Histone H3 (Phospho-T11) antibody (orb334714)

Description

Rabbit polyclonal antibody to Histone H3 (Phospho-



Western
blot
analysis of
Histone H3
(Pho...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T11 of human Histone H3 protein. The exact sequence is proprietary.

Target

HIST1H3A; HIST1H3B; HIST1H3C; HIST1H3D;
HIST1H3E; HIST1H3F; HIST1H3G; HIST1H3H;
HIST1H3I; HIST1H3J; HIST2H3A; HIST2H3C;
HIST2H3D; H3F3A; H3F3B

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[P68431](#), [P68433](#), [Q6LED0](#), [P84244](#), [P84228](#),
[Q71DI3](#), [P84245](#), [P84243](#)

Entrez

[8352](#), [8356](#), [8354](#), [8968](#), [15078](#), [8357](#), [319152](#),
[8358](#), [8355](#), [291159](#), [100361558](#), [8350](#), [15077](#),
[8353](#), [8351](#)

Dilution Range

WB: 1:500-1:1000

Expiration Date

12 months from date of receipt.