


Product Datasheet

LDHD antibody (orb334764)

Description	Rabbit polyclonal antibody to LDHD
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein corresponding to human LDHD.
Target	LDHD
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q86WU2 , Q7TNG8
Entrez	197257 , 52815
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of LDHD expression...