


Product Datasheet

LexA DNA Binding Region antibody (orb334765)

Description	Rabbit polyclonal antibody to LexA DNA Binding Region
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein corresponding to LexA DNA Binding Region.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	WB: 1:10000-1:50000
Expiration Date	12 months from date of receipt.



Western blot analysis of LexA DNA Binding...