

Product Datasheet

Anti-Human Serum Albumin Antibody-HRP labled
[Out of stock] (orb334847)

Description

Rabbit polyclonal antibody to

labeled

Species/Host

Rabbit

Reactivity

Human

Conjugation

HRP

Tested

ELISA

Applications

Immunogen

Human Serum Albumin

Target

ALB

Preservatives

Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.

Form/Appearance

Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P02768

Entrez

213

Dilution Range

ELISA: Use at an assay dependent dilution.

Expiration Date

12 months from date of receipt.