

Product Datasheet

Ovalbumin antibody (FITC) (orb334857)

Description	Rabbit polyclonal antibody against Ovalbumin conjugated to FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ELISA
Immunogen	Full length Ovalbumin
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Form/Appearance	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.