

Product Datasheet

Mouse 4-1BB /TNFRSF9 Protein (orb334874)

Description Mouse 4-1BB, Fc Tag (41B-M5258) is expressed from human 293 cells (HEK293). It contains AA Val 24 - ...

Reactivity Mouse

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target 4-1BB

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 46.7 kDa. The protein migrates as 50-66 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence NP_001070976.1

Purity ● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Source Mouse 4-1BB, Fc Tag (orb334874) is expressed from human 293 cells (HEK293). It contains AA Val 24 - Leu 211 (Accession # Q8R037-1).

MW 46.7 kDa

Uniprot ID [Q8R037](#)

NCBI [NP_001070976.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Mouse 4-1BB protein