

Product Datasheet

Human B7-H5 / Gi24 / VISTA Protein (orb334878)

Description Human B7-H5, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Phe 33 - Ala 194 (Acc...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target B7-H5

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 44.8 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.

Protein Sequence AAH20568

Purity 95%


Source Human B7-H5, Fc Tag (orb334878) is expressed from human 293 cells (HEK293). It contains AA Phe 33 - Ala 194 (Accession # AAH20568).

MW 44.8 kDa

Uniprot ID [Q9H7M9](#)

NCBI [AAH20568](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human B7-H5 protein