


Product Datasheet

Human CD163 Protein (orb334885)

Description	Scavenger receptor cysteine-rich type 1 protein M130 (CD163) is also known as hemoglobin scavenger r...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD163
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 110.4 kDa. The protein migrates as 135-140 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Q86VB7
Purity	95%
Source	Human CD163, His Tag (orb334885) is expressed from human 293 cells (HEK293). It contains AA Ser 42 - Ser 1045 (Accession # Q86VB7).
MW	110.4 kDa
NCBI	Q86VB7
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD163
protein