

Product Datasheet

Influenza A (H5N1) Hemagglutinin (HA) Protein (orb334898)

Description Neuraminidase (NA) and hemagglutinin (HA) are major membrane glycoproteins found on the surface of i...

Reactivity Virus

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Hemagglutinin (HA)

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 60.1 kDa. The protein migrates as 67-86 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.


Protein Sequence Asp 17 - Gln 531

Purity 95%

Source Influenza A [A/Hong Kong/483/97 (H5N1)] HA, His Tag (orb334898) is expressed from human 293 cells (HEK293). It contains AA Asp 17 - Gln 531 (Accession # AAC32099.1).

MW 60.1 kDa

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of Viral
Hemagglutinin...