




Product Datasheet

Cellular Apoptosis Susceptibility antibody (orb33490)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CSE1L
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	110417
Uniprot ID	P55060
NCBI	NP_001243064.1, NP_001307.2
Dilution Range	IHC-P: 1:100, IF/ICC: 1:100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of paraffin...